

Available Now

The first **FDA 510K
CLEARED** On-site
Saliva Test.

Oratect®

Made in USA

US Patent: 6248598

FDA 510K CLEARED

- Two test strip configurations of 6 drugs
- Rapid collection and testing in one step
- Results in 5 minutes
- User-friendly and non-invasive
- Avoids cross gender observation
- Reduces sample adulteration or contamination
- Includes transport tube for confirmation testing

Drugs of Abuse

Methamphetamines (ME)	50 ng/mL
Marijuana (TH)	40 ng/mL
Cocaine (CO)	20 ng/mL
Amphetamines (AM)	50 ng/mL
Opiates (OP)	40 ng/mL
Phencyclidine (PC)	10 ng/mL











Oratect® Reference Guide

The Oratect® Oral Fluid Drug Screen Device is a simple one-step test for the detection of drugs of abuse in oral fluid.

Oratect® Procedures -- Refer to package insert for detailed instructions and technical information.

<p>1. Remove the blue cap by holding the sides and pulling gently. This will expose the collection pad. Make sure there is a blue line present in each window area indicating an unused device.</p> 	<p>2. Open mouth and gently rub the collection pad inside mouth against cheek in a circular motion several (approximately 15-20) times. Make sure to keep head level.</p> 	<p>3. Gently rub the collection pad against the opposite cheek in circular motion several (approximately 15-20) times.</p> 
<p>4. Gently rub the collection pad on top of the tongue several (approximately 15-20) times. Do not chew, suck, bite or bend the collection pad.</p>  <p>5. Rub the collection pad underneath the tongue several (approximately 15-20) times.</p>	<p>6. Place the collection pad underneath the tongue for approximately 30 seconds to collect saliva. Instruct the donor to hold the device in place with hand. Sufficient amount of saliva collected is indicated by the flow of the blue lines. Repeat steps 2-6 until blue lines flow.</p> 	<p>7. Remove from mouth as soon as blue lines flow at both of the test windows. Re-cap the device.</p>  <p>8. Lay the device on a flat surface and read results in 5 minutes after removing the device from mouth. See Package Insert.</p>

Interpreting Oratect® Test Results

Invalid Result

When **no** colored band appears in the **CONTROL (C) region**, the test is **invalid** even if there is a band in the test region. **Repeat** the test with a new device.



Example Interpretation:

ME: Invalid
TH: Invalid
CO: Invalid
AM: Invalid
OP: Invalid
PC: Invalid

Negative Result

For each test, two colored bands should be observed:

- **One in the CONTROL (C) region**
- **One in the specific TEST region**

The color of the test band may be slightly darker or lighter than the control band. Any visible band that can be seen is a **negative** result.



Example Interpretation:

ME: Negative
TH: Negative
CO: Negative
AM: Negative
OP: Negative
PC: Negative

Presumptive Positive Result

A colored band at the **CONTROL (C) region** should be observed. When there is no colored band at the specific **TEST region**, the test is presumptive positive for that particular drug.



Example Interpretation:

ME: **Presumptive Positive**
TH: Negative
CO: Negative
AM: Negative
OP: Negative
PC: Negative