

SANI-CLOTH® HB

Germicidal Disposable Wipe



The industry's most diverse line of wipes meeting all hard, non-porous surface disinfecting needs.


Sani-Cloth® HB is an EPA registered hospital-grade disinfectant to accommodate the many situations and hard, non-porous environmental surfaces found in healthcare settings.

Sani-Cloth® HB kills bacteria, viruses, and fungi including:

Hepatitis B • HIV-1 • Herpes Simplex Type 2 • Methicillin Resistant *Staphylococcus aureus* (MRSA) • Vancomycin Resistant *Enterococcus* (VRE) • Influenza A2-Asian Virus • *Escherichia coli* strain O157:H7 *Pseudomonas aeruginosa* • *Salmonella choleraesuis* • *Staphylococcus aureus* • *Klebsiella pneumoniae* • Vaccinia Virus • *Candida Albicans* Adenovirus Type 2 • Trichophyton Mentagrophytes

Alcohol Free* Formula

- Bactericidal, fungicidal, and virucidal
- A one-step disinfectant that works in the presence of 5% organic soil
- Cleans and disinfects
- Tested effective against over 100 microorganisms**
- First germicidal wipe proven effective against the Human (not duck) Hepatitis B virus
- Contact time is 10 minutes or less**











SANI-CLOTH® HB Features and Benefits

- Unique Deep Well Lid – Provides ample space to store and access next wipe
- Natural Latex free, bleach free and phenolic free
- Pre-saturated with the correct amount of germicidal solution – No mixing, no measuring – saves valuable time
- Cleans and disinfects – All in one product
- Two (2) year shelf life from date of manufacture on product
- EPA registered; meets CDC and OSHA guidelines – Maintain compliance to protocols
- Large and extra-large non-folded wipes in canisters; and folded wipes in individual packets – More packaging and wipe size options for optimal placement and surface coverage
- Diamond-embossed material – is thick and strong to help clean soiled surfaces
- Sani-Bracket® 3-in-1 Wall and Mobile Equipment Attachment System – So versatile, place Sani-Cloth® HB germicidal disposable wipes on walls, IV poles, gurneys, anywhere they are needed most



SANI-CLOTH® HB

Germicidal Disposable Wipes

Description	Reorder No.	Wipe Size	Case Pack	Case Weight	Case Dimension L x W x H	Case Cube	Pallet TI/HI
WIPES							
 Sani-Cloth® HB Large	Q08472	6 x 6.75 in 15.25 x 17.14 cm	160/can 12 can/cs	24.75 lbs 11.22 kg	15.56 x 10.68 x 16 in 39.52 x 27.12 x 40.64 cm	1.54 cu ft 0.043 cu m	10/3
 Sani-Cloth® HB Extra Large	Q85484	8 x 14 in 20.32 x 35.56 cm	65/can 6 can/cs	14.3 lbs 6.48 kg	15.68 x 10.87 x 10.25 in 39.8 x 27.60 x 26.03 cm	1.01 cu ft 0.028 cu m	10/4
 Sani-Cloth® HB Individual Large	H21082	5 x 8 in 12.7 x 20.32 cm	50/box 10 box/cs	7.72 lbs 3.50 kg	13.12 x 7 x 9.56 in 33.32 x 17.78 x 24.28 cm	0.51 cu ft 0.014 cu m	19/5
 Sani-Cloth® HB Individual Extra Large	U22495	11.5 x 11.75 in 29.21 x 29.84 cm	50/box 3 box/cs	7.9 lbs 3.58 kg	13.93 x 4.81 x 10.5 in 35.38 x 12.21 x 26.67 cm	0.41 cu ft 0.011 cu m	24/3
SANI-BRACKET							
 Sani-Bracket® 3-in-1 (Sani-Cloth®) Labeled: Not for Use on Skin	P010900	<u>Use with</u> Q08472 Q85484	10 /cs	19.5 lbs	29.75 x 12.75 x 9.25 in 75.5 x 32.5 x 23.5 cm	2.03 cu ft	4/4
 Sani-Bracket® 3-in-1 (Sani-Cloth®) Labeled: Not for Use on Skin	P010901	Q08472 Q85484	1 /cs	1.75 lbs	6 x 6 x 8.25 in 15.25 x 15.25 x 21 cm	0.17 cu ft	42/4
 Sani-Bracket® 3-in-1 Pole Accessory Kit	P010600	P010900 P010901	1 /cs	1.65 lbs	5.5 x 5.25 x 3 in 14 x 13.5 x 7.5 cm	0.065 cu ft	72/8
 Sani-Bracket® for X-Large Box (Sani-Cloth)	P44210	U22495	1 /cs	0.60 lbs	5.25 x 5.25 x 5.25 in 13.5 x 13.5 x 13.5 cm	0.07 cu ft	63/9

*Contains less than 0.15% alcohol.

**Refer to the Technical Data Bulletin for complete list of microorganisms and contact times.

Professional Disposables International, Inc.
Two Nice-Pak Park
Orangeburg, NY 10962-1376
ISO 9001:2008 Certified
Fax: (845) 398-5347
1-800-999-6423 in USA

www.pdipdi.com

©2010 Professional Disposables International, Inc.
Sani-Bracket® and Sani-Cloth®
are registered trademarks of
Professional Disposables International, Inc.

